GE

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Contre-indications
Hydrofilm can be used for protection against bacterial penetration, external contaminants, and mechanical irritation of dry or slightly exuding, primary healing wounds, for almost healed epithelising wounds as well as for the securing of wound dressings, catheters and cannulas. Additionally, it can be used as a sealing element in negative pressure wound therapy.

Indications
Hydrofilm can be used for protection against bacterial penetration, external contaminants, and mechanical irritation of dry or slightly exuding, primary healing wounds, for almost healed epithelising wounds as well as for the securing of wound dressings, catheters and cannulas. Additionally, it can be used as a sealing element in negative pressure wound therapy.

Contraindications
Hydrofilm must not be used as a primary dressing on clinically infected, bleeding or highly exudative wounds.

Mode of application (see pictograms I-IV)
Clean the application area thoroughly and ensure that the skin is dry and free from creams and ointments. Hairly skin should be shaved before applying Hydrofilm. (I) Open the peel pack and take out the dressing. Hydrofilm is protected by 4 covers, numbered 1 to 4 in the order of application. Slightly remove the protective covers 1 and 2, then (II) place the Hydrofilm adhesive side on the wound. Fully remove the protective covers 1 and 2, gently press the entire film onto the skin, conforming it to the body contours. (III) Remove the protective cover 3. (IV) then remove the protective cover 4 and gently press Hydrofilm onto the skin. Hydrofilm can be removed almost painlessly by overstretching the material. To do this, loosen an edge and carefully take off the dressing by overstretching it.

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Side effects
Adverse reactions may occur in patients who are allergic or hypersensitive to acrylic adhesives.

Incident reporting
For a patient/user/third party in the European Union and in countries with identical regulatory regime (Regulation (EU) 2017/745 on Medical Devices), if, during the use of this device or as a result of its use, a serious incident occurs, please report it to the manufacturer and/or your authorized representative and to your national authority.

Product disposal
In order to minimize the risk of potential infection hazards, or environmental pollution, disposable components of the medical device should follow disposal procedures according to applicable and local laws, rules, regulations and infection prevention standards. Dispose of the medical device with household waste.

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Hydrofilm can be used for protection against bacterial penetration, external contaminants, and mechanical irritation of dry or slightly exuding, primary healing wounds, for almost healed epithelising wounds as well as for the securing of wound dressings, catheters and cannulas. Additionally, it can be used as a sealing element in negative pressure wound therapy.

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Mode of application (see pictograms I-IV)
Clean the application area thoroughly and ensure that the skin is dry and free from creams and ointments. Hairly skin should be shaved before applying Hydrofilm. (I) Open the peel pack and take out the dressing. Hydrofilm is protected by 4 covers, numbered 1 to 4 in the order of application. Slightly remove the protective covers 1 and 2, then (II) place the Hydrofilm adhesive side on the wound. Fully remove the protective covers 1 and 2, gently press the entire film onto the skin, conforming it to the body contours. (III) Remove the protective cover 3. (IV) then remove the protective cover 4 and gently press Hydrofilm onto the skin. Hydrofilm can be removed almost painlessly by overstretching the material. To do this, loosen an edge and carefully take off the dressing by overstretching it