

# Simplicity. Down to a science.

Compression and Fixation Bandage Systems



# Simple solutions for complex clinical challenges.

With our aging population, patients have at least one chronic condition and may have additional co-morbidities that can impact their quality of life. Healthcare providers across all care settings are stretched thin and costs keep mounting. Clinicians are charged with the mission to provide faster and more efficient care at lower costs and higher quality. Your primary objective is to link the patient's needs to clinical treatment goals. And when products are part of the treatment plan, you want simple solutions that are safe, comfortable and easy to use.

## **In this challenging environment, a simple formulary can make all the difference.**

HARTMANN puts the patient first. Our high-quality products are designed to meet the changing needs of patients as they transition from higher-acuity settings to more independent care. By offering the right combination of therapeutic products, we help take complexity out of the equation so you can focus on what is most important – your patient.

By choosing HARTMANN, you have a partner dedicated to helping you improve the patient experience. For nearly 200 years, we have been innovating and actively seeking to improve healthcare. Your dedicated team of product specialists will help you create a formulary with a range of products to achieve your clinical treatment goals, then provide onsite product training for your staff. Our solutions are designed to generate repeatable, consistent practice to help promote improved clinical outcomes – all priced to offer you superior value.

## **Patient-centric, high-quality products backed by a highly trained team. It's really that simple.**

# Protect. Treat.

The science of sim



## Solutions that **PROTECT**

Our products protect the skin to help reduce the risk for skin injury often related to adhesives and may be used to shield IV sites, drains and tubes from abrupt or unplanned removal.



# Support.

Support creates the ideal solution for every clinical need.



## Solutions that **TREAT**

We help you manage the process of healing for disease states such as vascular insufficiencies, lymphedema and post-operative edema.



## Solutions that **SUPPORT**

Our compression and fixation solutions support and help promote healing for joint and bone injuries, secure splints and function as positioning aides during surgical and interventional procedures.



Depend on our bandages to deliver accurate stretch, compression, comfort and support for optimal healing.

See our complete family of elastic, cohesive and multi-layer compression bandages, designed to:  
Protect. Treat. Support.





# Solutions that **Protect.**

Our products protect delicate skin and help prevent damage caused by tape or other extrinsic factors.



## Peha-haft® LF

- Unique blend of soft, skin-friendly gauze and cohesive fibers
- Ideal for delicate or traumatized skin
- Conforms exceptionally well in difficult-to-dress areas

**COHESIVE**



## Co-Wrap® LF

- Soft, breathable and absorbent cohesive bandage that adheres to itself without pins, clips or tape
- For use as a secondary bandage, for over-wrapping primary dressings and securing IV tubes
- Provides light, consistent and secure compression

**COHESIVE**



## Co-Lastic® LF

- Water resistant and lightweight
- Non-slip cohesive bandage that adheres to itself, not the skin
- Quick find starting edge and hand tears easily

**COHESIVE**



## Comperm® LF and Comperm®-PD LF

- Cotton and polyester knit bandage with a tubular construction;
- Provides 360 degrees of uniform, consistent support
- Commonly used for edema reduction, post-cast support and protect sensitive skin from adhesives or trauma

**TUBULAR BANDAGE**

# Solutions that Treat.

Premium product solutions that help in the treatment of vascular insufficiencies, lymphedema and post-operative edema.



## Co-Lastic® LF

- Water resistant and lightweight
- Non-slip cohesive bandage that adheres to itself, not the skin
- Quick find starting edge and hand tears easily

COHESIVE



## Deluxe® LF and Deluxe® 480® LF

- Durable orthopedic grade woven bandage with high cotton content
- Provides superior, consistent elasticity, compression and support\*
- Ideal for limb support following injury or surgery

ELASTIC



## EZe-Band® LF

- Knit weave and high cotton content feel provides controlled, consistent compression
- Provides superior absorbency and breathability\*
- Easy to apply and secure

ELASTIC



## Econo-Paste® and Econo-Paste® Plus Calamine

- Non-slip cotton stretch gauze with zinc oxide paste (also available with calamine)
- Stretches and conforms easily
- Vacuum-sealed foil pouch for longer shelf life and innovative plastic core for easy application

ZINC PASTE BANDAGE



## LoPress®

- Short stretch, inelastic cotton bandage
- Absorbent, breathable and allows flexibility of movement
- Durable enough to withstand repeated use and washing

INELASTIC



## TwoPress® 2, ThreePress® and FourPress®

- TwoPress® 2 includes a short stretch cotton padding layer and a cohesive wrap for consistent compression, easy application and improved compliance
- ThreePress® includes a cotton padding layer, high stretch layer with stretch indicator and a cohesive wrap
- FourPress® includes a cotton padding layer, short stretch layer, high stretch layer and a cohesive wrap

MULTI-LAYER COMPRESSION

\*In comparison to other HARTMANN USA products.

# Solutions that Support.

Our compression solutions support and help healing for joint and bone injuries, secure splints, use as positioning aides and support during surgical and interventional procedures.



## Deluxe® LF and Deluxe® 480® LF

- Durable orthopedic grade woven bandage with high cotton content
- Provides superior, consistent elasticity, compression and support\*
- Ideal for limb support following injury or surgery

ELASTIC



## EZe-Band® LF

- Knit weave and high cotton content feel provides controlled, consistent compression
- Provides superior absorbency and breathability\*
- Easy to apply and secure

ELASTIC



## Shur-Band® LF and Shur-Band® Ice Wrap LF

- Knit bandage for general support, compression or securing and treating varicosities, sprains and torn ligaments
- Self-closure flap fastens anywhere on the bandage for quick, easy and secure use without pins, clips or tape
- Ice Wrap designed with sewn pouch to hold standard size hot and cold packs

ELASTIC



## Co-Lastic® LF

- Water resistant and lightweight
- Non-slip cohesive bandage that adheres to itself, not the skin
- Quick find starting edge and hand tears easily

COHESIVE



## Medi-Rip® LF

- Convenient, multi-purpose self-adherent elastic support bandage
- Woven with 99% cotton yarn, yet tears easily without scissors
- Absorbent, lightweight and conforms easily

COHESIVE



## Esmark LF

- Provides maximum and consistent compression to extremities for exsanguination prior to surgery
- Provides consistent 400% elongation for ideal compression and ease of application
- For use in surgical procedures especially joint replacement

ELASTIC

# Ordering information

Order no.	Product name	Size	Pcs/ Pack	Pcks/ Case
932452	Peha-haft® LF	1" x 4.5 yds.		8
932441	Peha-haft® LF	1.5" x 4.5 yds.		90
932442	Peha-haft® LF	2.25" x 4.5 yds.		90
932443	Peha-haft® LF	3" x 4.5 yds.		120
932444	Peha-haft® LF	4" x 4.5 yds.		135
45120000	Co-Lastic® LF	1" x 5 yds.		30
45150000	Co-Lastic® LF	1.5" x 5 yds.		48
45200000	Co-Lastic® LF	2" x 5 yds.		36
45300000	Co-Lastic® LF	3" x 5 yds.		24
45400000	Co-Lastic® LF	4" x 5 yds.		18
45600000	Co-Lastic® LF	6" x 5 yds.		12
45510000	Co-Lastic® LF	1" x 5 yds. blue		30
45520000	Co-Lastic® LF	2" x 5 yds. blue		36
45110000	Co-Lastic® LF	1.5" x 5 yds. color pack		48
45210000	Co-Lastic® LF	2" x 5 yds. color pack		36
45310000	Co-Lastic® LF	3" x 5 yds. color pack		24
29100000	Co-Wrap® LF	1" x 5 yds.	24	4
29150000	Co-Wrap® LF	1.5" x 5 yds.		125
29200000	Co-Wrap® LF	2" x 5 yds.	12	8
29300000	Co-Wrap® LF	3" x 5 yds.	12	8
29400000	Co-Wrap® LF	4" x 5 yds.	12	8
25100000	Medi-Rip® LF	1" x 5 yds.	24	4
25200000	Medi-Rip® LF	2" x 5 yds.	12	8
25300000	Medi-Rip® LF	3" x 5 yds.	12	8
25400000	Medi-Rip® LF	4" x 5 yds.	12	8
25600000	Medi-Rip® LF	6" x 5 yds.	12	4
70000599	Medi-Rip® LF	1.5" x 5 yds.		250
25320000	Medi-Rip® LF	2" x 5 yds.	12	8
25330000	Medi-Rip® LF	3" x 5 yds.	12	8
25810000	Medi-Rip® LF	1.5" x 5 yds.		32
25820000	Medi-Rip® LF	2" x 5 yds.		24
25830000	Medi-Rip® LF	3" x 5 yds.		16
83010000	Comperm® LF	A, 1.63" x 11 yds.*		1
83020000	Comperm® LF	B, 2.5" x 11 yds.*		1
83030000	Comperm® LF	C, 2.75" x 11 yds.*		1
83040000	Comperm® LF	D, 3" x 11 yds.*		1
83050000	Comperm® LF	E, 3.5" x 11 yds.*		1

Order no.	Product name	Size	Pcs/ Pack	Pcks/ Case
83060000	Comperm® LF	F, 4" x 11 yds.*		1
83070000	Comperm® LF	G, 5" x 11 yds.*		1
83080000	Comperm® LF	J, 7" x 11 yds.*		1
83090000	Comperm® LF	K, 8.5" x 11 yds.*		1
83100000	Comperm® LF	L, 13" x 11 yds.*		1
47300000	Econo-Paste®	3" x 10 yds.	1	12
47400000	Econo-Paste®	4" x 10 yds.	1	12
47310000	Econo-Paste® + Cal	3" x 10 yds.	1	12
47410000	Econo-Paste® + Cal	4" x 10 yds.	1	12
38200000	Deluxe® LF	2" x 5.5 yds.	10	6
38300000	Deluxe® LF	3" x 5.5 yds.	10	6
38400000	Deluxe® LF	4" x 5.5 yds.	10	6
38600000	Deluxe® LF	6" x 5.5 yds.	10	6
38410000	Deluxe® LF	4" x 11 yds. (double-length)		6
38610000	Deluxe® LF	6" x 11 yds. (double-length)		6
27200000	Deluxe® LF Sterile	2" x 5.5 yds.	1	20
27300000	Deluxe® LF Sterile	3" x 5.5 yds.	1	20
27400000	Deluxe® LF Sterile	4" x 5.5 yds.	1	20
27600000	Deluxe® LF Sterile	6" x 5.5 yds.	1	20
28400000	Deluxe® LF Sterile	4" x 11 yds.	1	20
28600000	Deluxe® LF Sterile	6" x 11 yds.	1	20
39200000	Deluxe® 480® LF	2" x 5 yds.	10	6
39300000	Deluxe® 480® LF	3" x 5 yds.	10	6
39400000	Deluxe® 480® LF	4" x 5 yds.	10	6
39600000	Deluxe® 480® LF	6" x 5 yds.	10	6
59120000	EZe-Band® LF	2" x 5 yds.	10	6
59130000	EZe-Band® LF	3" x 5 yds.	10	6
59140000	EZe-Band® LF	4" x 5.5 yds.	10	6
59160000	EZe-Band® LF	6" x 5.5 yds.	10	6
59180000	EZe-Band® LF	4" x 11 yds.	12	3
59190000	EZe-Band® LF	6" x 11 yds.	12	2
59720000	EZe-Band® LF Sterile	2" x 5 yds.	1	36
59730000	EZe-Band® LF Sterile	3" x 5 yds.	1	36
59740000	EZe-Band® LF Sterile	4" x 5.5 yds.	1	36
59760000	EZe-Band® LF Sterile	6" x 5.5 yds.	1	36
59790000	EZe-Band® LF Sterile	6" x 11 yds.	1	24

\*Unstretched

## Ordering information - continued

Order no.	Product name	Size	Pcs/ Pack	Pcks/ Case
59520000	Shur-Band® LF	2" x 5 yds.	10	6
59530000	Shur-Band® LF	3" x 5 yds.	10	6
59540000	Shur-Band® LF	4" x 5 yds.	10	6
59560000	Shur-Band® LF	6" x 5 yds.	10	6
59570000	Shur-Band® LF (double-length)	4" x 10 yds.		6
59580000	Shur-Band® LF (double-length)	6" x 10 yds.		6
59920000	Shur-Band® LF Sterile	2" x 5 yds.	1	10
59930000	Shur-Band® LF Sterile	3" x 5 yds.	1	10
59940000	Shur-Band® LF Sterile	4" x 5 yds.	1	10
59960000	Shur-Band® LF Sterile	6" x 5 yds.	1	10
59880000	Shur-Band® LF Ice	4" x 30"		12
59890000	Shur-Band® LF Ice	6" x 30"		12
18310000	Esmark LF Sterile	3" x 9 ft.*		20
18320000	Esmark LF Sterile	3" x 12 ft.*		20
18410000	Esmark LF Sterile	4" x 9 ft.*		20
18420000	Esmark LF Sterile	4" x 12 ft.*		20
18610000	Esmark LF Sterile	6" x 9 ft.*		20
18620000	Esmark LF Sterile	6" x 12 ft.*		20

\*Unstretched

Order no.	Product name	Size	Pcs/ Pack	Pcks/ Case
43400000	FourPress®			
	Kit contains 4 bandages and 3 tape strips			8
	Padding bandage	4" x 3.8 yds.*		
	Crepe bandage	4" x 4.9 yds.		
	Compression bndg	4" x 9.5 yds.		
	Cohesive bandage	4" x 6.5 yds.		
43420000	ThreePress®			
	Kit contains 3 bandages and 3 tape strips			8
	Padding bandage	4" x 4.3 yds.*		
	Compression bndg	4" x 9.6 yds.		
	Cohesive bandage	4" x 6.5 yds.		
931687	TwoPress® 2			
	Kit contains 2 bandages			8
	Padding bandage	4" x 7.1 yds.*		
	Cohesive bandage	4" x 8.9 yds.		
42200000	LoPress®	2.3" x 5.4 yds.	1	12
42300000	LoPress®	3.1" x 5.4 yds.	1	12
42400000	LoPress®	3.9" x 5.4 yds.	1	12
42500000	LoPress®	4.7" x 5.4 yds.	1	12
42410000	LoPress® (double-length)	3.9" x 10.9 yds.	1	6

Contact your local HARTMANN representative, or visit us at  
[www.hartmannusa.com](http://www.hartmannusa.com) for more information.

HARTMANN USA, Inc.  
 Rock Hill, SC 29730  
 1-800-243-2294



Going further  
for health