

# Advanced Wound Management






## HydroTac® and HydroTac® Comfort

Hydroactive foam dressing

- Accelerates the healing process
- Stimulate epithelialization process due to growth factor accumulation
- Optimizes pH values for epithelialization
- Reduces the potential for sticking to wounds and supports painless dressing changes
- High retention under pressure, including when used in combination with compression
- Protects wound edges against exudate, reducing the potential for maceration

### Shapes:

HydroTac®  
Non-Adhesive



HydroTac®  
Comfort  
Adhesive




## PermaFoam® and PermaFoam® Comfort

High vapor permeable foam dressings




- Large pores and high absorption capacity effectively support wound cleansing
- Specific surface structure stimulates the formation of healthy granulation tissue
- Fast vertical absorption reduces the potential of maceration of wound perimeter
- High fluid handling capacity may reduce the need for dressing changes
- Can be used under compression because its high retention rate

### Shapes:

PermaFoam®  
Non-Adhesive




PermaFoam®  
Comfort  
Adhesive




### Indications:

- Treatment of light to medium exuding wounds during the granulation and epithelialization phases
- Treatment of dry or exuding wounds



### Compatibility with Compression Treatment.



### Indications:

- Management of chronic draining wounds including pressure ulcers, leg ulcers and diabetic foot ulcers
- For use on medium to heavy exuding wounds
- Ideal for use under compression bandaging systems



**CollaSorb®**  
Collagen dressing

- Composed of 90% native collagen and 10% calcium alginate
- Sterile, non-cytotoxic bovine-derived collagen
- Absorbs exudate up to 30 times its weight
- Easy to cut and apply



**Sorbalgon® Ag**  
Silver calcium alginate dressing

- Silver ions protect dressing from micro-organisms for up to 7 days
- Holds together wet or dry
- Forms gel to aid in autolytic debridement
- Highly absorbent

**Shapes:**

- 2" x 2"
- 4" x 4"
- 6" x 6"
- 1" x 12" Ribbon



**Sorbalgon®**  
Calcium alginate dressing

- Maintains contact with wound base reducing the potential of a moist chamber
- Avoids pain during dressing changes
- Protects wound edges against maceration
- Supports the healing of bleeding wounds
- May be left in place several days
- Maintains a moist wound environment

**Shapes:**

- 2" x 2"
- 4" x 4"
- 4" x 8"


**Shapes:**

- 2" x 2"
- 4" x 4"

**Indications:**


- Managing acute and chronic wounds with light to heavy drainage
- Heavily draining wounds, a foam may be used as a secondary dressing
- Moderately draining wounds, cover with a hydrocolloid or adhesive dressing
- Minimal exudate, cover with a non-woven retention tape

Silver elution of different silver dressings




Source: Data on file.

Tensile strength N/cm



Source: Data on file.


Absorbency ml/100 cm²



**Indications:**

- Ideal for treating chronic or acute moderate to heavily draining wounds

Absorbency ml/100 cm²



Source: Data on file.

Absorption in velocity of blood in seconds



Source: Terrill P., Sussman G. and Bailey M. Absorption of blood by moist wound healing dressings. Primary Intention 2003; 11(1):6-7

**Indications:**

- Ideal for treating chronic or acute moderate to heavily draining wounds



## Hydrofilm® and Hydrofilm® Plus

Self-adhesive film dressing

- Skin friendly
- Easy application while wearing gloves
- Rounded corners reduce potential roll-up and extend wear time
- Highly breathable film reduces maceration risk
- Extremely thin 25-micron film provides high-wearing comfort
- Absorbent pad does not stick to the wound (Hydrofilm Plus)
- Waterproof and germproof
- Sterile

Shapes:

HydroFilm®	HydroFilm® Plus
2.4" x 2.75"	2" x 2.8"
4" x 5"	3.5" x 4"
4" x 6"	3.5" x 6"
6" x 8"	
8" x 12"	

## FlexiCol®

Hydrocolloid dressing



- Forms a gel to prevent sticking to wound
- Thin beveled edges reduce potential roll-off
- High and quick fluid uptake
- Conforms to difficult-to-dress body contours
- Comfortable to wear
- Quick and easy to apply
- Secure initial adhesion

## Cosmopor®

Adhesive wound dressing



- Soft and gentle to the skin – can be used for patients with sensitive skin
- Conforms well to body contours and remains in place
- Rounded corners help reduce the need to change the dressing
- Good absorbency and cushioning effect
- Fluid repellent backing

## Omnifix®

Non-woven retention/abrasion tape



- Skin-friendly adhesive
- Especially soft and conformable
- Highly permeable to air and water vapor
- Helps reduce potential for maceration
- Provides reliable adherence
- Self-adhesive, nonwoven

## Peha-Haft®


Conforming cohesive bandage



- Conforming elastic bandage that sticks to itself and not the skin, hair or clothes
- Very kind to the skin
- Does not contain natural rubber latex, chlorine or acrylate
- No need for additional fixation
- High absorption and air permeable

## MVTR (Upright)

The higher the number the better the absorption.



Source: Data on file.


## Precaution:

- Should not be used as a primary dressing on clinically infected, bleeding or heavily secreting wounds

## Indications:

- Serves as a secondary dressing
- Protects against physical damage for dry primary wounds and nearly healed epithelial wounds
- Secures catheters or cannulae
- Reduces the risks for secondary infection
- Used for sterile post-operative care of slightly exuding wounds


## Shapes:



## Indications:

- For use in managing light to moderately draining chronic and acute wounds including: pressure ulcers, leg ulcers, second degree burns and superficial burns and abrasions

## Shapes:



## Indications:

- Protection of post-surgical wounds
- Wounds with low to moderate level of exudate
- Protection of wounds following wound cleaning

## Sizes:

- 2" x 2 yards
- 4" x 2 yards
- 2" x 10 yards
- 4" x 10 yards
- 6" x 10 yards



## Indications:

- Dressing retention
- Ideal for use on joints, round and angular parts of the body
- Securement of catheters, measuring instruments, probes and similar devices

## Sizes:

- 1" x 4.5 yards
- 1.5" x 4.5 yards
- 2.25" x 4.5 yards
- 3" x 4.5 yards
- 4" x 4.5 yards



## Indications:

- Safe dressing retention for use on joints and conical or rounded parts of the body
- Fixation of cannulae, padding material, etc.

# Ordering information

Order no.	Product name	Size	Pcs/ Pack	Pcks/ Case
900800	Cosmopor® Sterile	2" x 2.8"	50	20
900806	Cosmopor® Sterile	3.2" x 4"	25	22
900808	Cosmopor® Sterile	3.2" x 6"	25	14
900812	Cosmopor® Sterile	4" x 8"	25	12
900814	Cosmopor® Sterile	4" x 10"	25	8
900816	Cosmopor® Sterile	4" x 14	25	8
900820	Cosmopor® Sterile	4" x 4"	25	8
900823	Cosmopor® Sterile	6" x 6"	25	8
48600000	FlexiCol®	2" x 2"	20	1
48610000	FlexiCol®	4" x 4"	10	1
48620000	FlexiCol®	6" x 7" Sacral	5	1
48630000	FlexiCol®	3.5" x 4.75" Concave	10	1
48640000	FlexiCol®	4" x 4" Thin	10	1
48660000	FlexiCol®	6" x 6"	5	1
49200000	Sorbalgon®	2" x 2"	10	1
49210000	Sorbalgon®	4" x 4"	10	1
49230000	Sorbalgon®	4" x 8"	5	1
999609	Sorbalgon® Ag	1" x 12"	5	1
999610	Sorbalgon® Ag	6" x 6"	5	1
999611	Sorbalgon® Ag	4" x 4"	10	1
999612	Sorbalgon® Ag	2" x 2"	10	1
409427	PermaFoam®	2.5" Round	10	1
409401	PermaFoam®	4" x 4"	10	1
409402	PermaFoam®	4" x 5"	10	1
409405	PermaFoam®	6" x 6"	5	1
409406	PermaFoam®	8" x 8"	3	1
409425	PermaFoam®	4" x 4" Cavity	3	1
409426	PermaFoam®	3.2" x 3.2" Trach	10	1
409408	PermaFoam® Comfort	4.3" x 4.3"	10	1
409412	PermaFoam® Comfort	6" x 6"	5	1
409413	PermaFoam® Comfort	8" x 8"	3	1
409422	PermaFoam® Comfort	7" x 7" Sacral	3	1
409424	PermaFoam® Comfort	6.5" x 7" Concave	3	1
685755	Hydrofilm®	2.4" x 2.75"	10	1
685756	Hydrofilm®	2.4" x 2.75"	100	1
685757	Hydrofilm®	4" x 5"	10	1

Order no.	Product name	Size	Pcs/ Pack	Pcks/ Case
685758	Hydrofilm®	4" x 5"	100	1
685759	Hydrofilm®	4" x 6"	10	1
685760	Hydrofilm®	4" x 6"	50	1
685761	Hydrofilm®	6" x 8"	10	1
685762	Hydrofilm®	6" x 8"	50	1
685765	Hydrofilm®	8" x 12"	10	1
685770	Hydrofilm® Plus	2" x 2.8"	5	1
685771	Hydrofilm® Plus	2" x 2.8"	50	1
685772	Hydrofilm® Plus	3.5" x 4"	5	1
685773	Hydrofilm® Plus	3.5" x 4"	50	1
685774	Hydrofilm® Plus	3.5" x 6"	5	1
685775	Hydrofilm® Plus	3.5" x 6"	25	1
685849	HydroTac®	2.5" Round	10	1
685832	HydroTac®	4" x 4"	10	1
685833	HydroTac®	4" x 8"	3	1
685839	HydroTac®	6" x 6"	3	1
685843	HydroTac®	6" x 8"	10	1
685844	HydroTac®	8" x 8"	3	1
685810	HydroTac® Comfort	3" x 3"	10	1
685815	HydroTac® Comfort	5" x 5"	10	1
685817	HydroTac® Comfort	6" x 6"	3	1
685821	HydroTac® Comfort	6" x 8"	10	1
685822	HydroTac® Comfort	8" x 8"	3	1
89850000	CollaSorb®	2" x 2"	30	1
89860000	CollaSorb®	4" x 4"	30	1
932452	Peha-haft®	1" x 4.5 yds.		8
932441	Peha-haft®	1.5" x 4.5 yds.		90
932442	Peha-haft®	2.25" x 4.5 yds.		90
932443	Peha-haft®	3" x 4.5 yds.		120
932444	Peha-haft®	4" x 4.5 yds.		135
36200000	Omnifix®	2" x 2 yds. unstretched		24
36400000	Omnifix®	4" x 2 yds. unstretched		24
36210000	Omnifix®	2" x 10 yds. unstretched		24
36410000	Omnifix®	4" x 10 yds. unstretched		12
36610000	Omnifix®	6" x 10 yds. unstretched		12

To learn more about our products or to place an order contact  
**HARTMANN at 1-800-243-2294 or visit [hartmannusa.com](http://hartmannusa.com).**

To learn more about our products or to place an order contact  
HARTMANN at 1-800-243-2294 or visit [hartmannusa.com](http://hartmannusa.com).

HARTMANN USA, Inc.  
Rock Hill, SC 29730  
1-800-243-2294

[hartmannusa.com](http://hartmannusa.com)

XLIT 2698 Rev. 2 (0117)



Going further  
for health